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Stand der Information: November 2023

Referenzen

- Schulz P. et al. Biochemical characterization, stability, and pathogen safety of a new fibrinogen concentrate (fibryga®). *Biologicals* 2018;52:72–77
- Fachinformation November 2023
- Göbel K. et al. The Coagulation Factors Fibrinogen, Thrombin, and Factor XII in Inflammatory Disorders – A Systematic Review. *Front Immunol* 2018;9:1731–1745
- Ross C. et al. Pharmacokinetics, clot strength and safety of a new fibrinogen concentrate: randomized comparison with active control in congenital fibrinogen deficiency. *Journal of Thrombosis and Haemostasis*, 2018;16:253–261
- Lissitchkov T. et al. Efficacy and safety of a new human fibrinogen concentrate in patients with congenital fibrinogen deficiency: an interim analysis of a Phase III trial. *Transfusion* 2018; 58(2):413–422
- Roy A. et al. Efficacy of fibrinogen concentrate in major abdominal surgery – a prospective, randomized, controlled study in cytoreductive surgery for pseudomyxoma peritonei. *Journal of Thrombosis and Haemostasis* 2019. <https://doi.org/10.1111/jth.14665>. Published online October 2019.
- Callum J. et al. Effect of Fibrinogen Concentrate vs Cryoprecipitate on Blood Component Transfusion after Cardiac Surgery: The FIBRES Randomized Clinical Trial. *JAMA*. doi:10.1001/jama.2019.17312. Published online October 21, 2019.